



## Wrist Orthosis with Moldable Insert

### Features:

- Moldable plastic insert envelops greater than 50% of the circumference of the wrist when molded to ensure a perfect fit.
- A custom fit is quick, easy, and hygienic to achieve without fumes or dust.
- Strong, light-weight, and remoldable.
- Volumetric adjustability - Compensates for swelling or muscle atrophy throughout the healing process, making it an ideal fit.
- Latex-free

### Indications:

- Carpal Tunnel Syndrome
- Wrist sprains or strains
- Post-Operative care



Moldable plastic Insert

RIGHT	LEFT	SIZE	MEASURE CIRCUMFERENCE		
			HAND	WRIST	FOREARM
3920-SM-R	3920-SM-L	SMALL	7" - 8"	5" - 6"	6" - 7½"
3920-MD-R	3920-MD-L	MEDIUM	8" - 9"	6½" - 7½"	8" - 9½"
3920-LG-R	3920-LG-L	LARGE	9" - 10"	7½" - 8½"	8½" - 10"
3920-XL-R	3920-XL-L	X-LARGE	9½" - 10½"	8½" - 9½"	10½" - 12½"