



THERMO-MED SMART SUSPENSION BELT FOR FEMORAL PROSTHESIS

Features:

Neoprene thigh support with elastic distal straps made in conical shape adapted to the different circumferences of the thigh.

Function:

- Textured grip to secure your prosthesis.

Indications:

- Transfemoral prosthesis.
- Slight instabilities.
- Long-lasting painful processes.
- Stabilization during sports activities or situation of great effort.



LEFT	RIGHT	SIZE	CIRCUMFERENCE "A" IN CM	CIRCUMFERENCE "B" IN CM
5502-SM	5501-SM	SMALL / 2	85 - 100	40 - 50
5502-MD	5501-MD	MEDIUM / 3	85 - 110	48 - 58
5502-LG	5501-LG	LARGE / 4	105 - 125	56 - 66


