



3-Tex® 7104 KNEE BRACE WITH POLYCENTRIC JOINTS



Features:

Fabric knee support made from three-layer breathable elastic material. Composed of three layers – elastic microfibre fabric, polyurethane foam and terry cloth – which facilitates rapid dispersion of sweat, greater elasticity and a better fit. Features silicone bubble grip to prevent slippage. Open patella with silicone patella pad, metal hinges with polycentric joints, femoral and tibial straps for firmer fastening and calf opening. Latex-free.

Indications:

- Joint pain.
- Osteoarthritis and arthritis.
- Tendinitis.
- Mild mediolateral instability.
- Capsular ligament injuries.



ITEM	SIZE	SIZE	KNEE CIRCUMFERENCE IN CM
7104-XS	X-SMALL	1	32 - 34
7104-SM	SMALL	2	34 - 36
7104-MD	MEDIUM	3	36 - 38
7104-LG	LARGE	4	38 - 41
7104-XL	X-LARGE	5	41 - 44
7104-XXL	XX-LARGE	6	44 - 49

