



WRAP AROUND KNEE SUPPORT WITH MEDIOLATERAL SPIRAL STABILISERS

Features:

Fabric knee support in three-layer breathable elastic material. Made from elastic velour (knee area) and elastic fabric (strap area). Composed of three layers – elastic microfibre fabric, polyurethane foam and terry cloth – which facilitates rapid dispersion of sweat, greater elasticity and a better fit, and features silicone bubble grip to prevent slippage.

Open patella with patella pad, mediolateral straps and opening in the popliteal area.

Function:

Compression, support and stability.

Indications:

- Stabilization during sports activities or straining situations.
- Patients whose special morphology requires quick and easy placement.



ITEM	SIZE
7119-UNI	UNIVERSAL

