

9104 RODISIL CLOSED PATELLA KNEE BRACE

Closed elastic knee support with silicone pad and lateral stabilizers

Features:

Developed using three-dimensional woven fabric that provides greater transversal elasticity at the same time as it exerts the necessary compression in each area. This new woven fabric concept provides better hold on the skin and stability during movement and does not require any fastening system. Incorporated in it is a viscoelastic cushion in ring form which stabilizes the patella. It has flexible stabilizers on both sides which prevents it from slipping and helps to provide greater medial and lateral control. **The Rodisil provides a therapeutic compression of 20.33-36 mmHg.**

Function:

Support, compression, and adjustment. The silicone pad produces an intermittent micromassage that favours re-absorption of edemas and hematomas with every knee movement.


Indications:

- Tendomyopathy
- First or second degree ligament sprains
- Mild Osteoarthritis
- Post-traumatic and post-operative use
- Sensation of instability (patellar and/or medial-lateral)



Bilateral design fits both left and right knee

ITEM	SIZE	THIGH CIRCUMFERENCE IN CM (A)	CALF CIRCUMFERENCE IN CM (B)
9104-XS	X-SMALL	38 - 41	28 - 31
9104-SM	SMALL	41 - 44	31 - 34
9104-MD	MEDIUM	44 - 47	34 - 37
9104-LG	LARGE	47 - 50	37 - 40
9104-XL	X-LARGE	50 - 53	40 - 43
9104-XXL	XX-LARGE	53 - 56	43 - 46
9104-XXXL	XXX-LARGE	56 - 59	46 - 49
9104-4XL	4X-LARGE	59 - 62	49 - 52



A-Point of measurement: outline halfway up thigh
B-Point of measurement: outline halfway up calf
Height: 32 cm