

Wrist Lacer™

wrist support



8" wrist support

With the ease of application of the Med Spec Wrist Lacer, you don't have to worry about your patients being compliant. In addition, they will love the exceptional comfort and support the Wrist Lacer provides.

Ease of Application

The lace closure mechanism provides unparalleled ease of application. You only need one free hand to tension and fasten the lace closure. The stockinette sleeve prevents the thumb from catching the laces and provides padding between the laces and the patient's skin.

Adjustability & Immobilization

Malleable palmar and dorsal aluminum stays may be bent to the desired degree of flexion or extension. The blue palmar pad provides enhanced support by matching the contour of the depression in the palm.

Optimum Comfort

Made with perforated suede and a polypropylene felt liner, the Wrist Lacer allows moisture to evaporate away from the skin. The anatomically contoured design and the adjustable thumb strap provide a comfortable fit and allow unrestricted motion of the fingers.

Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Post operative care



10 1/2" wrist & forearm support

Size	8" black left	8" black right	8" beige left	8" beige right
x-small	223921	223931	223221	223231
small	223922	223932	223222	223232
medium	223924	223934	223224	223234
large	223925	223935	223225	223235
x-large	223926	223936	223226	223236

Size	10 1/2" black left	10 1/2" black right	10 1/2" beige left	10 1/2" beige right
small	223952	223962	223552	223562
medium	223954	223964	223554	223564
large	223955	223965	223555	223565
x-large	223956	223966	223556	223566



rev 7/07

Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com www.medspec.com suggested HCPCS: L3908

© 2007 Medical Specialties, Inc.