

Compressor™

wrist support



The compressor strap on this wrist support makes it an excellent choice for providing compression along the carpal region. In addition, it works in conjunction with the plastic stays to create a functional wrist support that limits the range of flexion and extension of the wrist. The patented strap retention system on the end of the compressor strap allows for ease of application

Coolflex™ material

breathable, four-way stretch, cushioning

Molded, high density plastic palmar stay

anatomical fit: resists flexion of wrist

Low density, plastic dorsal stay

resists extension of wrist

Compressor strap

patented strap retention system allows for ease of application



| Left | Right | Size | Wrist Circ. |
|--------|--------|--------|-------------|
| 223802 | 223812 | small | 5" - 6" |
| 223804 | 223814 | medium | 6" - 7" |
| 223805 | 223815 | large | 7" - 8" |