

GelFlex® Wrist



The GelFlex Wrist significantly reduces scarring at the incision site following carpal tunnel surgery, while permitting significant mobility of the hand and wrist. The contoured polymer gel pad is specially formulated to aid in healing and improving the skin's external forces, thereby reducing discomfort caused by hypersensitivity.

Features and benefits include:

A contoured polymer gel pad that is formulated with a medical grade mineral oil which gradually dissipates onto the skin to help reduce scar tissue following carpal tunnel surgery. The gel is clinically proven to be effective in improving the skin's elasticity while minimizing scarring. The soft gel pad (4½" long) also protects the surgical site from friction, pressure and vibration; ideal for return to work activities.

The Coolflex™ material does not contain any rigid stays and it is extremely comfortable, cool and low profile which promotes patient compliance.



Item #	Size	Wrist Circ.
223822 (left)	small	5" - 6"
223824	medium	6" - 7"
223825	large	7" - 8"
223826	x-large	8" - 9"
223832 (right)	small	5" - 6"
223834	medium	6" - 7"
223835	large	7" - 8"
223836	x-large	8" - 9"