

Motion Manager

Excellent choice for return to activity following wrist sprain or carpal tunnel.

Features & Benefits:

- * **Semi-flexible polyethylene stay:** controls wrist motion in functional range
- * **Stays are removable:** restriction in specific ranges can either be increased or decreased depending on individual patient needs
- * **Flip tab:** allows product to be worn on either right or left hand
- * **Neoprene base:** excellent cushioning and support
- * **Ventilated adjustable strap:** allows adjustment of compression around the carpals



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel: 1-866-334-4441
Fax: 1-888-484-4431

Item #	Size	Measurements
223901	X-Small	up to 5 3/4"
223902	Small	5 3/4" - 6 1/2"
223904	Medium	6 1/2" - 7 1/4"
223905	Large	7 1/4" - 8"
223906	X-Large	over 8"



Professional Orthopedic Products Providing Innovative Products of Unsurpassed Quality