

Ryno Lacer[®]

wrist & thumb support



Long Ryno Lacer

The Ryno Lacer provides superior immobilization of the CMC and MP joints of the thumb, yet it permits limited motion of the IP joint to allow for contact with the opposing fingers. The Ryno Lacer offers your patients the support they need and the comfort they demand.

Lace closure design

Allows for easy application

Radial, palmar and dorsal stays

Restrain wrist and thumb motion. Aluminum stays may be bent to the desired degree of flexion or extension of the wrist and abduction or adduction of the thumb.

Polypropylene liner and perforated suede

Allow moisture to evaporate away from the skin.

Comfort and compliance

The anatomically contoured design and the adjustable thumb strap provide a comfortable fit and allow unrestricted motion of the fingers.

Indications:

- DeQuervain's Tendonitis
- Gamekeeper's Thumb
- Post Cast Healing
- Selected Scaphoid Injuries



Short Ryno Lacer

Long Black Left	Long Black Right	Long Beige Left	Long Beige Right	Size
223981	223991	223781	223791	x-small
223982	223992	223782	223792	small
223984	223994	223784	223794	medium
223985	223995	223785	223795	large
223986	223996	223786	223796	x-large

Short Black Left	Short Black Right	Size
223881	223891	x-small
223882	223892	small
223884	223894	medium
223885	223895	large
223886	223896	x-large

