

EpiGel™

elbow support



A comfortable, low profile design that aids in the treatment of tennis elbow (lateral epicondylitis) or golfer's elbow (medial epicondylitis). The EpiGel elbow support provides relief with two adjustable gel pads that compress the extensor and flexor muscles without restricting circulation of the forearm. The soft gel pads also dampen vibration. The universal design fits a broad size range of patients.

item #	size	elbow circ.
223410	univ	8" - 14"

* measure forearm just below elbow

Features

- Two Gel pads
- Patented strap retention system
- Universal design
- Low profile design

Benefits

Compress extensor muscles to relieve lateral epicondylitis and flexor muscles to relieve medial epicondylitis without creating a tourniquet around the forearm. Gel pads also dampen vibration.

Allows for easy application.

Fits right or left arm. One size fits all.
Less likely to catch on any object.

EpiFoam™

elbow support

item #	size	elbow circ.
223420	univ	8" - 14"

* measure forearm just below elbow



Features

- Single foam pad
- Patented strap retention system
- Universal design
- Low profile design

Benefits

Compress extensor muscles to relieve lateral epicondylitis or flexor muscles to relieve medial epicondylitis.

Allows for easy application.

Fits right or left arm. One size fits all.
Less likely to catch on any object.

Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte NC 28227
 email: request@medspec.com www.medspec.com
 p: 800-334-4143 p: 704-694-2434 f: 800-694-9060 f: 704-694-9060
 suggested HCPC code: L3701

EpiGel is a trademark of Medical Specialties, Inc.

CE
rev0505