

POST-OP ELBOW BRACES



The Vantage Post-Op Elbow Braces

- Select From A Dual or A Single Upright System
- Two effective alternatives to casting in some of your elbow rehabilitation patients.
- Effective range of motion control utilizing Vantage's self-contained ROM hinge.
- These Vantage Elbow Braces were designed to enhance your elbow rehab management.
- Range of motion is controlled in 10 degree increments.
- Flexion - 120 to 0 degrees. Extension - 90 to 0 degrees.



Single Upright

Right Sm - 032051 Md - 032052 Lg - 032053

Left Sm - 032054 Md - 032055 Lg - 032056

Dual Upright Universal

Sm - 038SML Md - 038MED Lg - 038LNG

Sizing - Measuring from the axis of the elbow to the axis of the wrist



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel:1-866-334-4441
Fax:1-888-484-4431