

Med Spec® Shoulder Immobilizer



The Med Spec Shoulder Immobilizer is designed to provide your patients with the optimum level of comfort and shoulder stabilization. This device may be used to provide immobilization for Bankhart repairs, glenohumeral dislocations or subluxations, labral repairs, rotator cuff repairs, and soft tissue repairs or strains. The ease of application and the lightweight design of the Med Spec Shoulder Immobilizer improve patient compliance.

Features

- Soft padded shoulder strap
- Soft padded waist strap
- Bilateral
- Single plastic ring above hand

Benefits

- Comfortably distributes weight on healthy shoulder
- Restricts abduction and external rotation of the shoulder
- Fits right or left shoulder
- Allows for easier application of shoulder strap



item #	size	measurements*
183370	Universal	up to 15 1/2"
183371	x-small	up to 12 1/2"
183372	small	up to 13 1/2"
183374	medium	up to 14 1/2"
183375	large	up to 15 1/2"
183376	x-large	up to 17 1/2"

*Measure from Olecranon (tip of elbow) to MP joint (knuckles)



Medical Specialties, Inc.
4600-K Lebanon Rd.
Charlotte, NC 28227
p: 800.334.4143 f: 800.694.9060
p: 704.694.2434 f: 704.694.9060
e-mail: request@medspec.com
www.medspec.com
suggested HCPCS code: L3660