

Back-n-Black™



This lumbar-sacral support features two flexible plastic abdominal panels, a semi-rigid sacral pad, and an inelastic circumferential nylon strapping system that produces intra-abdominal pressure to reduce load on the lumbar spine. This patented design provides a surprising level of support and comfort.



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel:1-866-334-4441
Fax:1-888-484-4431

Item #	Size	Measurements
163201	x-small	22-26"
163202	small	27-31"
163204	medium	32-38"
163205	large	39-45"
163206	x-large	46-52"
163207	xx-large	53-59"
163208	xxx-large	over-59"

Professional Orthopedic Products Providing Innovative Products of Unsurpassed Quality