



Back-n-Black™ Low Profile

This lumbar support incorporates the same patented design features of the regular Back-n-Black. This support is ideal for short-statured or female patients as the vertical dimension is 2" shorter than the regular Back-n-Black. In addition, the soft plush lining allows this support to be worn comfortably next to the skin.



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel: 1-866-334-4441
Fax: 1-888-484-4431

Item #	Size	Measurements
163241	x-small	22-26"
163242	small	27-31"
163244	medium	32-38"
163245	large	39-45"
163246	x-large	46-52"
163247	xx-large	53-59"