

Med Spec® Sacroiliac (SI) Belt



With its unique design, the Med Spec Sacroiliac Belt helps stabilize and support the sacroiliac joints, thereby relieving stress and discomfort at these core weight bearing joints. The superior materials of the Med Spec Sacroiliac Belt allow for the brace to be worn comfortably underneath your clothing and to remain in its proper position throughout your daily activities. The Med Spec Sacroiliac Belt provides your patients with the stabilization they need and encourages their compliance with an optimum blend of comfort and support.



Features

- Inner latex-free elastic belt with plush lining
- Strategically placed SkinLoc™ pads
- Outer nylon belt
- Removable hi-density foam pad

Benefits

- provides compression and is comfortable enough to be worn directly next to the skin
- grip the skin or underlying clothing to prevent the belt from migrating
- does not stretch when moving or lifting so it can help stabilize the sacroiliac joints
- may be used to provide additional pressure on the sacrum

Item #	Size	Waist Dimen.
163302	small	26" - 31"
163304	medium	31" - 38"
163305	large	38" - 45"
163306	x-large	45" - 52"
163307	xx-large	52" - 59"



Sacrum Foam Pad



May be worn directly next to the skin



SkinLoc Pad



rev 4/06