

# Criss-Cross Lumbar Support

**CMO**  
CMO, Incorporated

Overlapping abdominal straps provide exceptional lumbar and abdominal support. The three separate closures provide a wide range of tension adjustment, assuring the best fit.

## **Easy To Put On**

Even for patients of below normal strength. Each of the three closures require little strength. Together, they add up to exceptional support.

## **Comfortable**

Criss-Cross conforms to the body with no stiff stays to poke and irritate. It is contoured in front for comfort and its elastic has a soft push on the skin side.

## **No Creep Or Crawl**

Criss-Cross is tailored to ensure even tension at top and bottom. Combined with the securing effects of its overlapping abdominal straps and optional heat moldable insert, annoying slippage is eliminated.



<b>WHITE</b>		
<b>Item #</b>	<b>Size</b>	<b>Measurements</b>
<b>0130 WHI XS/SM</b>	<b>x-small/small</b>	<b>24-34"</b>
<b>0130 WHI MD/LG</b>	<b>medium/large</b>	<b>35-45"</b>
<b>0130 WHI XL/2X</b>	<b>x-large/xx-large</b>	<b>46-60"</b>

<b>BEIGE</b>		
<b>Item #</b>	<b>Size</b>	<b>Measurements</b>
<b>0130 BEI XS</b>	<b>x-small</b>	<b>24-29"</b>
<b>0130 BEI SM</b>	<b>small</b>	<b>30-34"</b>
<b>0130 BEI MD</b>	<b>medium</b>	<b>35-39"</b>
<b>0130 BEI LG</b>	<b>large</b>	<b>40-45"</b>
<b>0130 BEI XL</b>	<b>x-large</b>	<b>46-52"</b>
<b>0130 BEI 2X</b>	<b>xx-large</b>	<b>53-60"</b>



**Professional Orthopedic Products**

5465 Queen Mary Rd  
Suite 430  
Montreal, Quebec,  
H3X 1V5  
Tel: 1-866-334-4441  
Fax: 1-888-484-4431

**Professional Orthopedic Products Providing Innovative Products of Unsurpassed Quality**