

Spina II Model 1200

SPINA® II Lightweight Boltaron® Plastic Spinal Orthosis An economical lightweight Boltaron® plastic chairback orthoses lined with a 3/8" loop compatible plush soft-foam. The apron front is the same 3/8" plush-foam material with vertical stays and durable 2" plastic loops. The 2" straps are made of double-sided Velcro® webbing, adjustable at both ends which provide a wide range of flexibility.



MALE

Item #	Size	Measurements
SPINAI 1200-MS	small	30-35"
SPINAI 1200-MM	medium	36-41"
SPINAI 1200-ML	large	42-47"
SPINAI 1200-MXL	x-large	48-51"

FEMALE

Item #	Size	Measurements
SPINAI 1200-FS	small	30-35"
SPINAI 1200-FM	medium	36-41"
SPINAI 1200-FL	large	42-47"
SPINAI 1200-FXL	x-large	48-51"



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel: 1-866-334-4441
Fax: 1-888-484-4431

Professional Orthopedic Products Providing Innovative Products of Unsurpassed Quality