

Spina II Model 1275

Spina® II Knight-Taylor Spinal Orthosis Economical lightweight Boltaron® plastic Knight-Taylor spinal orthoses lined with a 3/8" loop compatible plush soft-foam. The apron front is the same 3/8" plush-foam material with vertical stays and durable 2" plastic loops. The 2" straps are made of double-sided loop webbing, adjustable at both ends which provide a wide range of flexibility. Adjustable over-the-shoulder straps with foam axillary pads.



MALE

Item #	Size	Measurements
SPINAII 1275-MS	small	30-35"
SPINAII 1275-MM	medium	36-41"
SPINAII 1275-ML	large	42-47"
SPINAII 1275-MXL	x-large	48-51"

FEMALE

Item #	Size	Measurements
SPINAII 1275-FS	small	30-35"
SPINAII 1275-FM	medium	36-41"
SPINAII 1275-FL	large	42-47"
SPINAII 1275-FXL	x-large	48-51"



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel:1-866-334-4441
Fax:1-888-484-4431

Professional Orthopedic Products Providing Innovative Products of Unsurpassed Quality