

# 403J Series: Dynamic "J" Buttress Knee Support



Stabilization achieved via a dynamic tubular "J" buttress for superior patella support. Specify left or right. Adjustable support straps provide a customizable level of compression and increased stability. Additional support includes spiral stays and support straps above and below the knee.

\*\*\*Measure at mid-patella with leg at full extension.

Please Specify Right or Left leg.



Professional Orthopedic Products

5465 Queen Mary Rd  
Suite 430  
Montreal, Quebec,  
H3X 1V5  
Tel: 1-866-334-4441  
Fax: 1-888-484-4431

Item #	Size	Measurements
<b>403J-S</b>	<b>small</b>	<b>13"-14"</b>
<b>403J-M</b>	<b>medium</b>	<b>14"-15"</b>
<b>403J-L</b>	<b>large</b>	<b>15"-17"</b>
<b>403J-XL</b>	<b>x-large</b>	<b>16"-18"</b>
<b>403J-XXL</b>	<b>xx-large</b>	<b>18"-20"</b>

