

403 Series: Universal Knee Support



Stabilization achieved via a universal "horseshoe" buttress for medial, lateral, superior or inferior patella support. Additional support includes spiral stays and support straps above and below the knee.



***Measure at mid-patella with leg at full extension.



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel: 1-866-334-4441
Fax: 1-888-484-4431

Item #	Size	Measurements
403-S	small	13"-14"
403-M	medium	14"-15"
403-L	large	15"-17"
403-XL	x-large	16"-18"
403-XXL	xx-large	18"-20"
403-XXXL	xxx-large	20"-22"