

# 9958 Series: Patella Trac



**BODY GLOVE**  
PROTECT THE CORE



Active patellar subluxation reduction. Adjustable support straps provide medial or lateral support to limit patellar instability. Dynamic force maintains contact and stability throughout range of motion. McConnell taping replacement using dynamic/elastic straps to correct tracking problems. Low profile design allows for comfortable all day wear. The patella buttress applies consistent force throughout range of motion.



**Professional Orthopedic Products**

5465 Queen Mary Rd  
Suite 430  
Montreal, Quebec,  
H3X 1V5  
Tel: 1-866-334-4441  
Fax: 1-888-484-4431

Item #	Size	Measurements
<b>9958-XS</b>	<b>x-small</b>	<b>11.5"-13"</b>
<b>9958-S</b>	<b>small</b>	<b>13"-14.5"</b>
<b>9958-M</b>	<b>medium</b>	<b>14.5"-15.5"</b>
<b>9958-L</b>	<b>large</b>	<b>15.5"-16.5"</b>
<b>9958-XL</b>	<b>x-large</b>	<b>16.5"-17.5"</b>