

H408 Series: Hinged Stabilizer Knee Support



Pull on design for a quick, easy application. It provides medial and lateral support to enhance knee stability. Hinges prevent hyperextension and provide stabilization. A low profile design allows for comfortable all day wear. Open politeal eliminates bunching and irritation behind the knee.



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel:1-866-334-4441
Fax:1-888-484-4431

Item #	Size	Measurements
H408-S	small	13"-14"
H408-M	medium	14"-15"
H408-L	large	15"-17"
H408-XL	x-large	16"-18"
H408-XXL	xx-large	18"-20"
H408-XXXL	xxx-large	20"-22"

