

ADVANTAGE II ACL READY FIT KNEE BRACE

Construction:

Lightweight, streamlined aluminum frame, sturdy to hold the fitting. The Advantage II will be better at keeping it's fit than other off-the-shelf ACL braces.

High tech hard powdercoating provides appearance and outstanding durability.

One of, if not the lowest profiles of any functional ACL brace on the market today, ready-fit or custom.

Comfort:

The Advantage II is ready to fit , easily adjusted to each of your patient's needs. Greater adjustability provides greater comfort... and performance!

BREATHEOPRENE **SEVERtec** cuff liners allow heat and moisture to dissipate. And straps are removable and hand washable. Important for your most active patients.

The Advantage II's low profile was designed with all day wear in mind, it can be worn all day... every day. Cuff geometry is designed to offer enhanced soft-tissue control and at the same time comfort.

Cost-Control:

Designed for the new world of health care... high quality at a reasonable price. Vantage has utilized the latest manufacturing processes to keep the quality high and the costs low.

For a quote on price you'll need to ask your Vantage representative or call us directly at 800-952-7223.

Control:

Advantage II's anatomically correct VectorGlide hinge mimics the natural motion of the human knee. Providing self-contained range of motion control setting.... no wedges, stops or other parts to lose. These settings allow for greater versatility throughout rehabilitation.

The Advantage II utilizes a rigid frame which provides a fulcrum for support. For your ACL patient the frame work with the VectorGlide hinge to reduce A-P translation and provide range of motion control.

For your MCL/LCL patient, adjustable condylar pads provide medial-lateral support.



Sizing Chart

X-Small	14 3/4" or less
Small	15" to 17 3/4"
Medium	18" to 20 1/2"
Large	20 3/4" to 22"
X-Large	22 1/4" to 23 1/2"
XX-Large	23 3/4" or more

Please Specify: Left or Right

6" Above Joint Line