

# The Vantage Full Leg Foam System



High grade HyperFoam wraps provide comfort and durability. Integrated circumferential straps make application a snap. Strong formable uprights give optimal contour and varus/valgus control.



Professional Orthopedic Products

5465 Queen Mary Rd  
Suite 430  
Montreal, Quebec,  
H3X 1V5  
Tel:1-866-334-4441  
Fax:1-888-484-4431

Item #	Size	
<b>C13000GREY</b>	<b>26"</b>	<b>Universal</b>
<b>C13000SGREY</b>	<b>23"</b>	<b>Universal</b>
<b>C13DELAGREY</b>	<b>20"</b>	<b>Universal</b>



*Reliable sturdy Range of Motion hinge that allows adjustment in 10 degree increments.*