

Universal Foam Knee Immobilizer



This lightweight foam knee immobilizer provides a wide range of fittings. A key benefit is the reduction of migration due to the five elastic circumferential straps holding it firmly in place. Malleable medial and lateral stays provide customization for additional support. The heavy posterior stay keeps the leg in full extension.



Professional Orthopedic Products

5465 Queen Mary Rd
Suite 430
Montreal, Quebec,
H3X 1V5
Tel: 1-866-334-4441
Fax: 1-888-484-4431

Item #	Size	Measurements
223660	16" universal	up to 24½"
223665	18" universal	up to 26½"
223670	20" universal	up to 29"
223675	22" universal	up to 31½"
223680	24" universal	up to 33¼"