

Gripper™

hinged knee brace



The Gripper™ Hinged Knee Brace aids in the treatment of lateral or medial collateral ligament sprains. The brace may also be used post-operatively or during sports to provide medial and lateral stability and to protect against hyperextension of the knee. The Gripper's unique design resists distal migration. The wraparound anterior closure allows for ease of application and an optimum fit.

The Gripper is available in two different base materials: Neoprene and CoolFlex™. The neoprene version provides warmth and superior compression. The CoolFlex version is breathable and hypoallergenic. Both versions of the Gripper are latex free.

Features

Internal SkinLoc™ patch

Outer closure straps

Steel hinge with reinforced pocket

Open popliteal

Bilateral design

Benefits

Works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby promoting patient compliance.

Straps are made of a quality elastic and woven nylon materials. The superior durability of the straps allows for numerous applications without fraying or losing closure strength. The high tension elastic is positioned over the quadriceps and calf muscles to allow for expansion during muscle contraction.

Hinged pocket is lined with ballistic nylon material to improve durability

Minimizes material bunching during flexion to improve comfort.

Fits either right or left knee.



Neoprene



CoolFlex

Neoprene	Size	Knee Circ.	CoolFlex
117011	x-small	10" - 12"	117031
117012	small	12" - 14"	117032
117014	medium	14" - 16"	117034
117015	large	16" - 18"	117035
117016	x-large	18" - 20"	117036
117017	xx-large	20" - 22"	117037
117018	xxx-large	over 22"	117038

Medical Specialties, Inc.
 4600-K Lebanon Rd.
 Charlotte, NC 28227
 p: 800.334.4143 f: 800.694.9060
 p: 704.694.2434 f: 704.694.9060
 e-mail: request@medspec.com
 www.medspec.com
 suggested HCPCS code: L1820
 © 2006 Medical Specialties, Inc.



*measure circumference of knee at mid patella