

Dynatrack™ patella stabilizer



An easy to apply design, the DynaTrack patella stabilizer aids in the treatment of patellofemoral dysfunction. Each size fits either the right or left knee. The internal buttress can be positioned to apply either lateral or medial pressure on the patella. The amount of pressure can be adjusted by changing the tension on the straps leading from the internal buttress and the tension on the outer wrap. A large popliteal opening permits high degrees of flexion with no discomfort.

The DynaTrack with neoprene provides warmth and superior compression. The DynaTrack with CoolFlex™ material is breathable and provides excellent compression.



Neoprene



tensioning of the
internal buttress



large popliteal opening



CoolFlex

| Neoprene | Size | Knee Circ. | CoolFlex |
|----------|-----------|------------|----------|
| 117001 | x-small | 10" - 12" | 117021 |
| 117002 | small | 12" - 14" | 117022 |
| 117004 | medium | 14" - 16" | 117024 |
| 117005 | large | 16" - 18" | 117025 |
| 117006 | x-large | 18" - 20" | 117026 |
| 117007 | xx-large | 20" - 22" | 117027 |
| 117008 | xxx-large | over 22" | |

*measure circumference of knee at mid patella

Medical Specialties, Inc.
4600-K Lebanon Rd.
Charlotte, NC 28227
p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com
www.medspec.com
suggested HCPCS: L1800 + L2795
U.S. Patent # 7,004,919

©2006 Medical Specialties, Inc.
Dynatrack is a trademark of Medical Specialties, Inc.



rev 3/07