



ASO[®] with Stays Ankle Stabilizer



The Choice of ChampionsSM

With the new improvements to the ASO with Stays, the ASO maintains its leadership position in the ankle brace market. Give your athletes and patients the patented and proven technology of the ASO with Stays. The leader continues to stay several steps ahead...

Plastic Stays

Removable plastic stays allow for enhanced proprioceptive support.

Stabilizing Straps

The stabilizing Straps are now angled to better contour above the malleoli for enhanced comfort and better migration resistance. Nylon stabilizing straps form a complete figure-eight to help stabilize the subtalar and talocrural joints.

Elastic Cuff

Circumferential woven elastic cuff keeps laces and stabilizing straps secure. In addition, it works in concert with the upper boot to resist the inversion or eversion forces.

CoolFlex Achilles Pad

CoolFlex[™] padding has been increased on the inside of the boot along the Achilles tendon to provide greater comfort.

Hidden Seam

There is no longer a seam under the foot. The seam is now hidden along the rear portion of the boot to improve comfort.



Indications:

- Inversion or Eversion Ankle Sprains
- Medial or Lateral Ankle Instability
- Subtalar Joint Instability

Suggested HCPCS:
L1902

Men's Shoe Size	6-7	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	xxl
Women's Shoe Size	5-6	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	m	l
Wide Foot	s	m	l	l	xl

Bilateral design fits left or right foot

Black	White	Size	Ankle Circum.*
264031	223181	x-small	10" - 11"
264032	223182	small	11" - 12"
264034	223184	medium	12" - 13"
264035	223185	large	13" - 14"
264036	223186	x-large	14" - 15"
264037	223187	xx-large	15" - 16"

*measure ankle circumference according to diagram below



Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
 e-mail: request@medspec.com www.medspec.com U.S. Patent # 5,067,486
 ASO is a registered trademark of Medical Specialties, Inc.
 © 2009 Medical Specialties, Inc.

