



ASO®

EVO™

Ankle Stabilizer



The next evolution in ankle bracing technology

The ASO EVO ankle stabilizer is an evolutionary step forward in ankle protection with its unique stabilizing straps and cuffs that incorporate the three main elements of an ankle taping application: stirrups, heel lock (figure-eight), and circumferential closure. This new design enables the ASO EVO to comfortably provide exceptional ankle support and stability in the treatment of ankle injuries and to reduce the severity and frequency of future ankle injuries.

Internal Plastic Cuff

Dynamic cuff system serves the dual purpose of (1) redirecting the stirrup strap to initiate a figure-eight around the ankle and (2) applying circumferential compression around the distal fibula and tibia to help stabilize the distal tibiofibular joint.

Stirrup Strap

Inner stirrup strap is lined with Skinloc™ material to better capture the calcaneus (heel) and resist inversion or eversion of the ankle.

Stabilizing Straps

Form a figure-eight by first crossing over the anterior region of the tibia, fibula, and talus to help stabilize the talocrural joint and then crossing under the calcaneus to help stabilize the subtalar joint.



Inner cuff lining removed in the above photo to reveal plastic cuff system.

Indications:

- Inversion or Eversion Ankle Sprains
- Syndesmotic Ankle Sprain (High Ankle Sprain)
- Medial or Lateral Ankle Instability
- Subtalar Joint Instability

Bilateral design fits left or right foot

Internal plastic cuff and stabilizing straps are superimposed in the above photo for illustrative purposes.

Item #	Size	Ankle Circum.*
264091	x-small	10" - 11"
264092	small	11" - 12"
264094	medium	12" - 13"
264095	large	13" - 14"
264096	x-large	14" - 15"
264097	xx-large	15" - 16"



*measure ankle circumference according to diagram above

Men's Shoe Size	6-7	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	xxl
Women's Shoe Size	5-6	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	m	l
Wide Foot	s	m	l	l	xl

Medical Specialties, Inc.
4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com www.medspec.com suggested HCPCS: L1906 patent pending

ASO is a registered trademark of Medical Specialties EVO is a trademark of Medical Specialties, Inc.

© 2008 Medical Specialties, Inc.



rev 1/08