

NEW

ASO® Flex-Hinge™

Ankle Stabilizer



The Choice of ChampionsSM

The ASO Flex-Hinge combines the proven technology of the ASO's patented strapping system with an articulated plastic shell for superior ankle stabilization.



Internal Plastic Shell

The upper portion of the internal plastic shell is incorporated into the eyelet lace closure to create an inner cuff around the tibia and fibula which resists separation of the tibia and fibula to help treat syndesmosis ankle sprains [high ankle sprains]

Internal plastic shell with hinge is built into ballistic nylon boot and flexes to conform to the anatomy of the ankle, which greatly enhances patient comfort and compliance

Indications:

- Inversion or eversion ankle sprains
- Syndesmosis ankle sprains
- Ankle/subtalar instability
- Medial/lateral ankle instability

Stabilizing Straps

Figure eight nylon stabilizing straps and elastic cuff work in conjunction with the internal plastic shell to capture the calcaneus and resist inversion or eversion motion

Bilateral design fits left or right foot

Internal plastic shell is superimposed in the above photo for illustrative purposes

Item #	Size	Ankle Circum.
264072	small	11" - 12"
264074	medium	12" - 13"
264075	large	13" - 14"
264076	x-large	14" - 15"
264077	xx-large	15" - 16"

*measure ankle circumference according to diagram below



Men's Shoe Size	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15	15+
Narrow Foot		s	m	l	xl
Regular Foot	s	m	l	xl	xxl
Wide Foot	m	l	xl	xxl	xxl
Women's Shoe Size	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14	
Narrow Foot		s	m	l	
Regular Foot	s	m	m	l	
Wide Foot	m	l	l	xl	