

# ASO MAX™



The ASO MAX is similar in design and material to the regular ASO, except for the addition of Skinloc™ material lining the inside of the brace. This material allows the brace to be worn with or without a sock. When the ASO MAX is worn directly in contact with the skin, it provides maximum support and excellent comfort.

## ASO® sizing chart

To determine the proper circumference of the ankle, measure the distance around the heel.



Professional Orthopedic Products

5465 Queen Mary Rd  
Suite 430  
Montreal, Quebec,  
H3X 1V5  
Tel:1-866-334-4441  
Fax:1-888-484-4431

Item #	Size	Measurements
<b>264041</b>	<b>x-small</b>	<b>10-11"</b>
<b>264042</b>	<b>small</b>	<b>11-12"</b>
<b>264044</b>	<b>medium</b>	<b>12-13"</b>
<b>264045</b>	<b>large</b>	<b>13-14"</b>
<b>264046</b>	<b>x-large</b>	<b>14-15"</b>
<b>264047</b>	<b>xx-large</b>	<b>15-16"</b>