

NEW

Phantom™ Dorsal Night Splint



Your patients with plantar fasciitis can sleep easy with the Phantom™ Dorsal Night Splint. It is so comfortable they will hardly realize they are wearing a night splint.

Dorsal Design

Comfortable low profile design enhances patient compliance and there are no hard surfaces to contact opposing leg during sleep. In addition to less bulk, the breathable CoolFlex™ material is cooler to wear.

Adjustable Aluminum Stays

Fully padded internal dorsal stays may be bent to increase or decrease dorsiflexion to the desired level.



Non-skid Sole

Non-skid material sewn to sole provides better traction when patients need to ambulate during the night.

Indications:

- Plantar Fasciitis
- Achilles Tendonitis

Elastic Tensioner

Elastic material sewn under midfoot provides additional tension on plantar fascia.



Item #	Size	Women's shoe size	Men's shoe size
264122	small	4 1/2 - 7	5 1/2 - 8
264124	medium	7 1/2 - 9 1/2	8 1/2 - 10 1/2
264125	large	10 - 12	11 - 12 1/2
264126	x-large	12 1/2 +	13 +

Bilateral design fits left or right foot

Medical Specialties, Inc.
4600-K Lebanon Rd.
Charlotte, NC 28227

p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060

e-mail: request@medspec.com www.medspec.com
suggested HCPCS: L4396 Patent pending

© 2006 Medical Specialties, Inc.

CE
rev 11/06