


# Restorative™ Hand

Effectively and comfortably treat any wrist/hand/finger lost range and tension from fresh stroke, brain or spinal cord injuries to moderate deformities and pain.



## Features:

- Flex design encourages joints to relax.
- Heat mold to provide just low load passive stretch on the joints, then remold to gradually correct joints to normal alignment.
- Provides prolonged low load passive stretch to relengthen shortened tissue.
- Kydex base extends to end of thumb.
- Comfortably supports osteo and rheumatoid arthritic joints and treats.
- Prevents deformities if treated in early stages.



Size	Width of MCP	Wrist to Fingers
X-Small	2.75" - 3"	Up to 5"
Small	Up to 3"	Up to 6.5"
Medium	Up to 3.5"	Up to 7"
Large	Up to 3.75"	Up to 8"



*Professional Orthopedic Products*

Toll Free Tel: 1-866-334-4441  
Toll Free Fax: 1-888-484-4431  
Email: [profortho@bellnet.ca](mailto:profortho@bellnet.ca)  
[www.poportho.com](http://www.poportho.com)