

NEW

Med Spec® Padded Knee Sleeve



When you want to provide your athletes and patients with the therapeutic benefits that neoprene offers, choose the **Med Spec Padded Knee Sleeve**.

Premium Grade Neoprene

Two sided nylon with 1/8" neoprene provides compression, warmth, and durability.

Anterior Pad

Full 1/8" thick anterior pad offers enhanced protection and comfort.

Open Patella

Open patella relieves patella pressure and assists in optimum positioning.

Closed Popliteal

Closed popliteal and anatomical design maintain uniform compression around the knee.

Part	Size	Knee Circ.*
117061	x-small	12" - 13"
117062	small	13" - 14"
117064	medium	14" - 15"
117065	large	15" - 16"
117066	x-large	16" - 18"
117067	xx-large	18" - 20"
117068	xxx-large	20" - 23"

*measure circumference of knee at mid patella

CE
rev 10/08



Medical Specialties, Inc.
4600-K Lebanon Rd. Charlotte, NC 28227
p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com www.medspec.com
suggested HCPCS: L1825
© 2008 Medical Specialties, Inc.