

AB100 BOXIA FOOT SPLINT PLUS



Characteristics:

- The **Boxia Plus** allow patients to be autonomous and independent because of the new semi-automatic, one-handed closing system so they can enjoy a better quality of life in their day-to-day activities.

One supramalleolar part made of micro-perforated breathable fabric that wraps around the ankle. It features a pad in the back that hugs the Achilles tendon, provides maximum comfort and distributes pressure. There are also relief points offering a massage effect on the soft parts. The central anchoring system can be adjusted to the ankle perimeter to guarantee the most adequate position for each patient.

- An elastic traction band keeps the ankle in a neutral position and connects the supramalleolar part to the “S” hooks, or the hook to the footwear.
- The anti-slip hook attached to the footwear and a strap that is secured to the laces prevents undesirable movement.
- Four pairs of “S” hooks for greater stability and varum-valgum foot control.



Benefits:

- Keeps the ankle in a functional position to prevent foot collapse;
- Enables a more natural and efficient walking pattern

Indications:

- Weakness or inability to dorsiflex the ankle (clubfoot) as a result of neurological injuries, congenital alterations or muscular disorders.
- Flaccid paralysis associated with hemiparesis.



ITEM	SIZE	PERIMETER IN CM
AB100-SM	Small / 1	14 - 20
AB100-MD	Medium / 2	20 - 26
AB100-LG	Large / 3	26 - 32