



## ADVANTAGE III WALKER

- Constructed of durable, malleable aluminum uprights that can be formed for a customized fit;
- Premium softgood liner provides additional comfort;
- Ultra low heel height and rocker sole facilitates a normal gait pattern;
- Superior strength, extremely lightweight fracture walker.



### Indications:

- Soft tissue injuries,;
- Grade 2 and 3 sprains;
- Stable fractures;
- Post-operative, trauma, or rehab use.



Size	High Top	Low Top	Men's Shoe Size	Women's Shoe Size
SMALL	T32022	T32122	4½ - 7	6 - 8
MEDIUM	T32024	T32124	7½ - 10	8½ - 11½
LARGE	T32026	T32126	10½ - 12½	12 - 13½
X-LARGE	T32028	-----	13+	14+