



## PEDIATRIC WALKING BOOT

### Features:

Designed specifically for children and teens, the Non-Air Pediatric Walking boot provides essential support and comfort for those recovering from ankle or foot injuries. It's lightweight polymer construction and thoughtful design features make it an ideal choice for young patients needing stability and mobility during the healing process.

### Function:

#### Lightweight Polymer Construction:

Crafted from premium injection molded material, this walking boot offers superior support while minimizing recovery strain.

#### Bilateral Design:

Featuring a bilateral design, this walking boot can be worn on either the left or right foot.

#### Rocker Sole:

The innovative rocker sole design enables a natural walking gait. This feature promotes a more comfortable and natural walking motion, reducing strain on the injured area. The thermoplastic rubber sole is equipped with traction grooves, ensuring stability and preventing slips and falls.

#### Comfortable Polyester Velvet Felt Interior:

The interior provides a plush and comfortable fit that reduces pressure points which assists with overall patient compliance.

### Indications:

- Soft tissue injuries;
- Sprains and strains;
- Immobilization of fractures;
- Post-Operative, trauma, or rehab use.



|                 | SMALL          | MEDIUM         | LARGE             |
|-----------------|----------------|----------------|-------------------|
| Product Code    | 12-100-PEDI-SM | 12-100-PEDI-MD | 12-100-PEDI-LG    |
| Age             | 1 - 2.5        | 2.5 - 5        | 5 - 8             |
| Foot Length     | Up to 6"       | Up to 6" - 7"  | Up to 7" - 8"     |
| Leg Length      | Up to 9"       | Up to 11"      | Up to 14"         |
| Child Shoe Size | Up to 8.5      | Up to 8.5 - 11 | Up to 11.5 - 12.5 |