

NEW AirShift™ OA Knee Brace



The AirShift OA Knee Brace uses the power of air to comfortably relieve pain for those patients with unicompartmental osteoarthritis. New lighter weight and lower profile design allows for full flexion of the knee and less bunching in the popliteal.

Air Bladder System

Inflatable air bladder allows patients and practitioners to easily adjust the level of unloading. Double layer air bladder consists of an inner vinyl layer and outer nylon layer for greater durability. High quality air valve prevents leakage. Repositionable condylar gel pad comfortably distributes the air bladder's unloading force on the knee. Air bladder system and single upright hinge are applied opposite the affected side of the knee.



Wraparound Design

Wraparound closure allows for ease of application and an optimum fit. Sizing options enable you to inventory fewer SKUs while still fitting a wide range of patient sizes.



Spiral Strapping System

Padded spiral straps work in concert with the air bladder to unload the osteoarthritic compartment of the knee using a 3-point leverage system. Distal spiral strap and inner Skinloc™ material capture the upper calf muscle to resist distal migration.

Indications:

- Unicompartmental osteoarthritis of the knee
- Postoperative unicompartmental articular and meniscal cartilage repairs

Suggested HCPCS: L1843 or L1851

Left Medial Comp. or Right Lateral Comp.	Right Medial Comp. or Left Lateral Comp.	Size	6" Above	Mid-Patella	6" Below
117184	117194	Small/Medium	15" - 19"	12" - 16"	12" - 16"
117186	117196	Large/X-Large	19" - 25"	16" - 20"	16" - 20"
117188	117198	2X/3X-Large	25" - 33"	20" - 24"	20" - 24"

Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
 e-mail: request@medspec.com www.medspec.com
 Patent pending. AirShift is a trademark of Medical Specialties, Inc.
 © 2015-2022 Medical Specialties, Inc.