

# AKS™ Knee Support

## with Metal Hinges & Straps



The AKS Knee Support offers your patients the unique ability to adjust compression and fit from top to bottom.



### Adjustable Knee Sleeve

The distal and proximal openings allow for easy application and ability to fit various patient shapes. In addition, compression can be adjusted to the desired level.

### Adjustable Patellar Buttress

Soft horseshoe buttress is repositionable and works in concert with central sleeve region to provide circumferential compression around the knee and focal compression around the patella.

### Medial and Lateral Metal Hinges & Straps

Metal hinges and outer closure straps provide medial and lateral stability. Internal Skinloc™ pad and metal hinges function synergistically to resist distal migration and vertical bunching. Condylar pads enhance comfort.

### Indications:

- Arthritis of the knee
- Bursitis of the knee
- Chondromalacia patella
- Lateral or medial collateral ligament sprains
- Lateral or medial instability of the knee
- Patellar tendinitis

Bilateral design fits left or right knee.

Neoprene version provides greater compression and warmth.

CoolFlex version is breathable and lightweight for superior comfort.

Neoprene	Size	Knee Circ.*	CoolFlex
117111	X-Small	10" - 12"	117121
117112	Small	12" - 14"	117122
117114	Medium	14" - 16"	117124
117115	Large	16" - 18"	117125
117116	X-Large	18" - 20"	117126
117117	XX-Large	20" - 22"	117127
117118	3X-Large	over 22"	117128

Suggested HCPCS: L1820



\*measure circumference of knee at mid patella

Medical Specialties, Inc.  
 4600-K Lebanon Rd.  
 Charlotte, NC 28227  
 p: 800-334-4143 f: 800-694-9060  
 p: 704-694-2434 f: 704-694-9060  
 e-mail: request@medspec.com  
 www.medspec.com  
 Patent pending  
 ©2010-2023 Medical Specialties, Inc.