

# AKS™ Knee Support with Plastic Hinges



The AKS Knee Support offers your patients the unique ability to adjust compression and fit from top to bottom.



### Adjustable Knee Sleeve

The distal and proximal openings allow for easy application and ability to fit various patient shapes. In addition, compression can be adjusted to the desired level.

### Adjustable Patellar Buttress

Soft horseshoe buttress is repositionable and works in concert with central sleeve region to provide circumferential compression around the knee and focal compression around the patella.

### Medial and Lateral Plastic Hinges

Plastic hinges enhance the brace's ability to conform to the knee during flexion. Internal Skinloc™ pad and plastic hinges function synergistically to resist distal migration and vertical bunching.

### Indications:

- Arthritis of the knee
- Bursitis of the knee
- Chondromalacia patella
- Patellar tendinitis



Bilateral design fits left or right knee.  
Neoprene version provides greater compression and warmth.  
CoolFlex version is breathable and lightweight for superior comfort.

Neoprene	Size	Knee Circ.*	CoolFlex
117091	X-Small	10" - 12"	117101
117092	Small	12" - 14"	117102
117094	Medium	14" - 16"	117104
117095	Large	16" - 18"	117105
117096	X-Large	18" - 20"	117106
117097	XX-Large	20" - 22"	117107
117098	3X-Large	over 22"	117108

Suggested HCPCS: L1810



\*measure circumference of knee at mid patella

Medical Specialties, Inc.  
4600-K Lebanon Rd.  
Charlotte, NC 28227  
p: 800-334-4143 f: 800-694-9060  
p: 704-694-2434 f: 704-694-9060  
e-mail: request@medspec.com  
www.medspec.com  
©2010-2023 Medical Specialties, Inc.