

Back-n-Black™

Back Support



The patented design features of the Back-n-Black™ provide your patients with the support and comfort they need in the treatment of lumbar spine related disorders.

Anterior and Posterior Panels

The semi-rigid plastic anterior panel and plastic & foam posterior panels work in concert with the inelastic nylon strapping system to increase the intra-abdominal pressure, thereby decreasing the load on the intervertebral discs and surrounding musculoskeletal structures of the lumbar spine.

Nylon Strapping System

The new larger handles on the elastic belt and nylon strapping system allow for ease of application and tightening to create the desired level of compression.

Breathable Elastic Panels

The breathable elastic panels on the interior offer adjustability in fitting various abdominal shapes and enhance patient compliance with their low profile design.

Indications:

- Chronic or acute low back pain
- Ligamentous sprains & muscle strains of the lumbar spine
- Postural support



Item #	Size	Waist Circumference
163201	x-small	22" - 26"
163202	small	27" - 31"
163204	medium	32" - 38"
163205	large	39" - 45"
163206	x-large	46" - 52"
163207	xx-large	53" - 59"
163208	xxx-large	over 59"

PDAC Verified HCPCS: L0626 or L0641

