

CMC-X™ Thumb Stabilizer



The CMC-X™ Thumb Stabilizer provides your patients with pain relief from carpometacarpal (CMC) joint arthritis while allowing them to continue the activities they enjoy.

“X” Strapping System

The “X” strap is easily adjusted to position the first metacarpal (base of thumb) in the desired degree of abduction to help stabilize the CMC joint while still allowing the thumb to have functional mobility.

Adjustable Stay Pod

The adjustable stay pod works in concert with the “X” strap to provide gentle pressure and help stabilize the CMC joint. The stay pod contains a malleable aluminum stay and is easily bent and repositioned by a medical provider to provide a customized fit.

CoolFlex™ Material

The CMC-X is primarily made from breathable, hypoallergenic CoolFlex™ material for comfort and compression. The CMC-X does not contain any natural latex rubber or neoprene.



Size	Left	Right	Measurements
x-small	224451	224461	under 5½”
small	224452	224462	5½” - 6¼”
medium	224454	224464	6¼” - 7”
large	224455	224465	7” - 7¾”
x-large	224456	224466	over 7¾”

*measure circumference of wrist at narrowest point

Indications:

- CMC joint osteoarthritis (OA)
- Medial subluxation of the CMC joint

Suggested HCPCS: L3923 or L3924



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143
 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
 e-mail: request@medspec.com www.medspec.com U.S. Patent # 10,842,661
 © 2016-2024 Medical Specialties Inc.

