



CP04 POST OPERATIVE SHOE



Features:

- Non-slip rocker sole;
- Velcro closure system accommodates deformities or bulky post-operative dressings;
- Insole can be replaced with active off-loading insole for diabetic ulcers;
- Optional toe cap is easily removed;
- Made from breathable materials

Indications:

- Post-operative care of the foot and toes, in cases where it is not yet possible for the patient to wear conventional footwear due to venous return problems (oedema).
- Diabetic foot, plantar ulcers and other foot disorders when the off-loading insole is inserted.



Bilateral design fits both left and right foot

ITEM	SIZE	European Footwear Size
CP04-XS	X-SMALL	33 - 35
CP04-SM	SMALL	36 - 38
CP04-MD	MEDIUM	39 - 40
CP04-LG	LARGE	41 - 42
CP04-XL	X-LARGE	43 - 45