

Cubital Tunnel Brace



The new design features a mesh sleeve for easier application and improved comfort while still offering a low-profile adjustable angle night splint that comfortably holds the elbow in a neutral position to reduce tension and pressure on the ulnar nerve.



Ventral Design

Low-profile design enhances patient compliance during sleep. No rigid materials contact the head of the ulnar or the olecranon process.



Adjustable Aluminum Stay

Positions the elbow in a neutral position, yet the angle may be adjusted by bending the aluminum stay.

V-Strap™ Retention System

Provides easier application by preventing the proximal and distal closure straps from disengaging with the plastic buckles. Middle strap may be removed if desired for various treatment protocols.

Indications:

- Cubital tunnel syndrome
- Ulnar nerve entrapment at the elbow

Suggested HCPCS: L3762

Item #	Size	Circumference*
223444	Small/Medium	6" - 15"
223446	Large/X-Large	10" - 19"

*measure circumference of upper arm at mid bicep

