



Dynamic hip abduction
orthosis for children

by Dr. J. Visser

Dynamic hip abduction orthosis for children by Dr. J. Visser

Intended use is to support maturing of the hip with hip dysplasia and used after re-positioned hip dislocation with infants. The orthosis fixes the hip joint in a 90° flexion and 60° abduction position..

Indications

- Hip dysplasia (ICD10: Q65.8)
- Hip joint luxation (ICD10: S73.00)
- Hip joint subluxation (ICD10: T14.3)

Product functions

- Limiting range of motion
- The thigh bones are positioned in a way to help mature and secure the femoral head into the acetabulum to prevent unwanted movements and thereby hip (Sub) Luxation

Features

- Hip joints; 90° Extension stop and 30° Abduction stop
- Sideways adjustable hip joints
- Foam padded Thigh Cuffs with Velcro closure
- Chest strap with Velcro closure
- Material: back part made of plastic with polyethylene foam lining
- Color: white



Size chart

Article number	Size
2990	universal