

## E-240 STRONG SHOULDER SUPPORT



### Features:

Consists of two cotton elastic straps with independent rings for each side of the shoulders. The straps exert a gradual pull on the shoulders aiding in the correction of Kyphotic posture.

### Indications:

- Dorsalgia.
- To correct poor posture.



ITEM	SIZE	BELOW CHEST IN CM
E-240-1	XS/1	75 - 85
E-240-2	SM/2	80 - 92
E-240-3	MD/3	85 - 97
E-240-4	LG/4	92 - 104
E-240-5	XL/5	98 - 110
E-240-6	XXL/6	104 - 120

