

2024 Product Catalogue



**POP**

***Professional Orthopedic Products***

***Providing Innovative Products  
of Unsurpassed Quality***

**Proudly Celebrating over 28 Years Serving Canadians**





*Professional Orthopedic Products*

Professional Orthopedic Products (POP) has been a leading provider of Orthopedic, Sports Medicine, and Restorative products throughout all of Canada since 1996.

We import and distribute high-quality brands from around the world such as Medical Specialties, Orliman, Allard/CSUS, Basko, Restorative Medical, Scotts/CMO, and Elite Orthopedics just to name a few.

Physicians, orthotists, physiotherapists, and athletic therapists recommend our products for their patients. Consequently, our products are used by professional athletes from the NHL, the NBA and the CFL as well as hospitals, universities, and the Canadian Armed Forces.

From lower extremity to upper extremity to such well-known products as the ASO Ankle Brace, at Professional Orthopedic Products, we are dedicated to providing the highest quality of products and exceptional customer service.

If you would like more information about our products or are interested in a free in-service at your location, please call our office toll-free at 1-866-334-4441 or e-mail us at [csr@poportho.com](mailto:csr@poportho.com). Our dedicated and knowledgeable representatives are at your service. You can also visit us online at [www.poportho.com](http://www.poportho.com).

Thank you for your interest in Professional Orthopedic Products and we hope to include you among our many satisfied customers.



Professional Orthopedic Products

## Table of Contents:

Ankle	1 - 2
Arm	8 - 9
Back	7
Calf	3
Elbow	8 - 9
Foot	2 - 3
Groin	7
Hand	9 - 12
Knee	4 - 6
Shoulders	8
Thigh	6
Thumb	13
Torso	7
Wrist	9 - 12

<b>Province</b>	<b>Representative</b>	<b>Contact Information</b>
Alberta/Saskatchewan	Jason Orme Sales Representative	(403) 803-0884 jason.orne1973@gmail.com
British Columbia	Stephen Bourne Territory Manager	(778) 228-0234 sbourne@poportho.com
Maritimes	Scott Hubley Sales Representative	(902) 471-0301 scott@poportho.com
Newfoundland	Patrick Downey Sales Representative	(709) 687-1653 patrpdw@aol.com
Ontario South-West	Darryl Nequest Sales Manager	(519) 241-2705 dnequest@poportho.com
Quebec	Dominic Fafard Sales Representative	(514) 791-4779 dominic@poportho.com
Quebec	Jonathan Trottier Sales Representative	(514) 880-7561 jonathan@poportho.com

## EVO® - Ankle Stabilizing Orthosis



Patented design enables the EVO® to comfortably provide greater support than the original ASO® and better protection against high ankle sprains. Features an internal stirrup, & external stabilizing straps and dynamic cuff.

The EVO® is bilateral, available in black or white, and is available in sizes ranging from X-Small to XX-Large

## EVO® Speed Lacer™ - Ankle Stabilizing Orthosis



Has the same internal stirrup, external stabilizing straps, and dynamic cuff as the EVO® with the addition of a speed lace closure that offers a faster application.

The EVO® Speed Lacer™ is bilateral, available in black or white, and is available in sizes ranging from X-Small to XX-Large

## ASO® - Ankle Stabilizing Orthosis with Stays



Has the same great features of the original ASO® along with the addition of removable plastic stays which allow for enhanced proprioceptive support.

Available in black or white and is available in sizes ranging from X-Small to XX-Large

## ASO® Speed Lacer™ - Ankle Stabilizing Orthosis



Has the same great features as the original ASO® with the addition of the speed lace closure mechanism that offers a faster application.

The ASO® Speed Lacer™ is bilateral, available in black or white, and is available in sizes ranging from Pediatric to XX-Large

## EVO® Quatro™ - Ankle Stabilizing Orthosis



Low profile version of the EVO®. Ideal for athletes or patients that do not tolerate bulky materials inside their shoes. Features an internal stirrup, & external stabilizing straps and dynamic cuff.

Available in sizes ranging from X-Small to XX-Large

## OS6240 Elastic Ankle Support with Gel Pads



Pull-on ankle sleeve made of a breathable stretch knit fabric. Includes two viscoelastic gel pads which redistribute compression from the malleolus to the surrounding tissue area. The massaging effect of the compressive fabric and padding stimulates faster reabsorption of bruises and swellings.

Available in sizes ranging from Small to Large

## TOB-500N Crossover Elastic Ankle Support



Sock-type ankle support made of soft elastic. Features figure 8 straps that allow adjustable compression and support.

Available in black or beige in sizes ranging from X-Small to X-Large

## TGO490 TheraGo Functional Elastic Support



The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. In addition, it optimizes blood circulation and facilitates body temperature regulation.

Available in sizes from X-Small to X-Large

## ACE902 Actius Elastic Ankle Support



Made from breathable elastic, it is lined with a highly resistant soft and stretchable fabric that provides greater comfort.

Indications include: Mild injuries and sprains; Dislocations; Support and protection during sporting and work-related activities.

Available in sizes from Small to Large

## Phantom Dorsal Night Splint



The comfortable low-profile design enhances your patient compliance. Padded internal dorsal stays to increase or decrease dorsiflexion. Non-skid material on sole provides better traction when ambulating.

Available in sizes from Small to XX-Large

## RMI SafeBoot II



Flanged heel elevates the heel and relieves the pressure thereby preventing development of heel ulcers and allows for healing of pre-existing ulcers. The SafeBoot II keeps patient's foot in dorsiflexion which prevents drop foot or heel cord contractures. The rotator bar, when extended either medially or laterally, inhibits any internal or external rotation at the hip and therefore prevents hip contractures.

Available in adult sizes Small to X-Large as well as Pediatric sizes 3" to 8".

## Elements Stretch-UP



Sock style treatment for plantar fasciitis. Keeps tension on the plantar fascia so it heals in a stretched position at night.

Bilateral design available in Sizes from Small to Large

## AB100 Boxia Plus



The Boxia Plus allow patients to be autonomous and independent because of the new semi-automatic, one-handed closing system so they can enjoy a better quality of life in their day-to-day activities. It features a pad in the back that hugs the Achilles tendon, provides maximum comfort and distributes pressure. There are also relief points offering a massage effect on the soft parts. The central anchoring system can be adjusted to the ankle perimeter to guarantee the most adequate position for each patient.

Sizes ranging from Small to Large.

## FSX Pneumatic Walker



Durable, semi-rigid shell helps support the limb while providing full-shell protection. The integrated inflation system for separate medial and lateral air bladders means no need for a separate air pump. The low-profile rocker sole enables a natural gait pattern. Features a covered toe pad for additional protection.

Available in low-top and high-top versions in sizes ranging from Small to Large.

## Advantage III - Fixed Walker



An extremely lightweight and superior strength fracture walker constructed of durable, malleable aluminum uprights that can be formed for a customized fit with a premium softgood liner that provides additional comfort. The ultra-low heel height and rocker sole facilitates a normal gait pattern. Available in both a high-top and low-top version.

Sizes range from Small to Large

## M2 Night Splint



Helps stretch the plantar fascia with the adjustable bilateral dorsiflexion straps.

Bilateral design available in sizes Small to Large

## TGO484 TheraGo Functional Calf Support



The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. In addition, it optimizes blood circulation and facilitates body temperature regulation.

Available in sizes from X-Small to X-Large

## OS6804 Elastic Calf Support



Manufactured in a breathable stretch knit fabric that is very strong and soft for added comfort; its seamless construction means no uncomfortable rubbing.

Indications include: Indicated for muscle disorders and tibial protection; Support and protection for athletic and work activities

Available in sizes from Small to Large

## Gripper™ - 12" and 16" Hinged Knee Brace



Provides medial and lateral stability and resists distal migration due to an internal Skinloc™ pad. Wrap-around anterior closure allows for ease of application and optimal fit. Open popliteal minimizes material bunching.

Available in CoolFlex and Neoprene in 12" or 16".

Available sizes range from X-Small to 5X-Large

## Patellavator™ - Knee Orthosis



Has a unique design that applies variable pressure to the patellar tendon without creating a tourniquet around the leg or irritating the popliteal. Low profile design eliminates irritation to the opposite leg.

Available in sizes Small to 2X-Large

## DynaTrack™ Plus - Patella Stabilizer



The DynaTrack™ Plus Patella Stabilizer contains medial and lateral plastic hinges which enhance the brace's ability to conform to the knee during flexion. Dynamic buttress, lined with our patented Skin-Loc™ material applies continuous pressure to the patella during full range of motion of the knee. Open Popliteal permits a high degree of flexion with no discomfort.

Available in both CoolFlex and Neoprene material.

Available in ranging sizes from X-Small to 3X-Large.

## DynaTrack™ with Hinges - Patella Stabilizer



The DynaTrack™ with Metal Hinges has the same patented design as the DynaTrack Plus with the addition of medial and lateral metal hinges for extraordinary stability.

The DynaTrack™ with Metal Hinges must be purchased in a left or right configuration. The left configuration applies pressure to the lateral aspect of the patella on the left knee and the right configuration applies pressure to the lateral aspect of the patella on the right knee.

Available in sizes from X-Small to 2X-Large

## AKS™ with Plastic Hinges - Knee Support



Distal and proximal openings allow for easy application to fit various patient shapes. Soft horseshoe buttress is repositionable. Internal SkinLoc™ pad and plastic hinges function synergistically to resist distal migration and vertical bunching.

Sizes X-Small to 3X-Large.

## AKS™ with Metal Hinges and Straps



Distal and proximal openings allow for easy application to fit various patient shapes. Soft horseshoe buttress is repositionable. Internal SkinLoc™ pad and metal hinges function synergistically to resist distal migration and vertical bunching.

Available in sizes ranging from X-Small to 3X-Large.



**9104 Rodisil Closed Patella Knee Brace**

Three-dimensional fabric provides a therapeutic compression of 20.33mm to 36mm Hg. Viscoelastic buttress helps stabilize as well as massage the patella to reduce swelling. Metal side stabilizers help to prevent distal migration.

Available in sizes from X-Small to 4X-Large

**9106 Rodisil Closed Patella with Metal Hinges**

Three-dimensional fabric provides a therapeutic compression of 20.33mm to 36mm Hg. Bilateral polycentric metal hinges control medial/lateral instabilities. Silicone pad massages the patella to reduce swelling.

Available in sizes from X-Small to 4X-Large

**3-Tek® 7104 Knee Brace with Polycentric Joints**

Fabric knee support made from three-layer breathable elastic material. Composed of three layers – elastic microfibre fabric, polyurethane foam and terry cloth – which facilitates rapid dispersion of sweat, greater elasticity and a better fit. Features silicone bubble grip to prevent slippage. Open patella with silicone patella pad, metal hinges with polycentric joints, femoral and tibial straps for firmer fastening and calf opening.

Sizes X-Small to 2X-Large.

**K80032 Post-Op Adjustable Knee Brace**

Designed to provide controlled range of motion for patients recovering from knee surgery or those who have knee injuries or instabilities. Comfortable, simple to apply and extremely lightweight. Telescoping bars are easily adjusted for a personalized fit. Drop lock features allows brace to be locked out with the push of a button. Easy to use ROM hinge with quick adjusting flexion and extension stops to control and limit knee flexion and extension.

Available in one Universal size.

**ACE803 ELASTIC KNEE SUPPORT w/ STABILIZERS**

Made from breathable elastic, it is lined with a highly resistant soft and stretchable fabric that provides greater comfort.

Indicated for Arthritis and osteoarthritis, mild bruising, inflammatory processes and support during sporting and work-related activities.

Available in sizes ranging from Small to Large

**ACE802 ELASTIC KNEE SUPPORT**

Made from breathable elastic, it is lined with a highly resistant soft and stretchable fabric that provides greater comfort.

Indicated for Arthritis and osteoarthritis, mild bruising, inflammatory processes and support during sporting and work-related activities.

Available in sizes ranging from Small to Large

## Orliman Elastic Knee with Lateral Stabilizers



Manufactured in a breathable stretch knit material. Includes a viscoelastic pad and lateral metal stabilizers for improved medial/lateral control. The compressive fabric and viscoelastic pad improves circulation and helps to reduce swelling.

Available in two lengths:  
OS6211 - 24 cm  
OS6212 - 28 cm

Sizes range from Small to Large

## MedSpec Padded Knee Sleeve



When you want to provide your athletes and patients with the therapeutic benefits that neoprene offers, choose the MedSpec Padded Knee Sleeve.

The two-sided nylon with 1/8" neoprene provides compression, warmth, and durability with a full 1/8" thick anterior pad for enhanced protection and comfort.

The open patella relieves patella pressure and assists in optimum positioning while the closed popliteal and anatomical design maintain uniform compression around the knee.

Available in sizes X-Small to 3X-Large.

## TGO480 Functional Elastic Knee Support with Strips



The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. In addition, it optimizes blood circulation and facilitates body temperature regulation.

Available in sizes from X-Small to X-Large

## TGO487 Functional Elastic Knee Support with Joints



The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. In addition, it optimizes blood circulation and facilitates body temperature regulation.

Available in sizes from X-Small to X-Large

## TGO481 TheraGo Functional Elastic Thigh Support



The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. In addition, it optimizes blood circulation and facilitates body temperature regulation.

Available in sizes from X-Small to X-Large

## OS6801 Elastic Thigh Support



Manufactured in a breathable stretch knit fabric that is very strong and soft for added comfort; its seamless construction means no uncomfortable rubbing.

Indications include: Prevention of tears of the quadriceps and abductor muscles; Prevention of tears of the ischiotibial and adductor muscles; Support and protection for athletic and work activities.

Available in sizes from Small to Large.

## Loving Comfort Maternity Support



Maximum comfort for women with moderate symptoms. The Loving Comfort Maternity Support is made from lightweight ventilated elastic that is cool and feather soft to the skin.

Available in pre-pregnancy dress sizes from Small to Large

## Reinforced Hernia Truss with Velcro Fastening



Elastic belt with Velcro fastening, anatomic cushions extended with straps under the buttocks and fixed to the sides.

Indicated for bilateral reducible inguinal hernias; compression and reduction of hernias; maximum containment and comfort.

Models available:

B-110:Bilateral / D-111:Right / I-112:Left

Sizes range from 70 cm to 120 cm of waist circumference.

## Back-n-Black™ Back Support



Patented design includes anterior and posterior panels combined with inelastic nylon strapping increases intra abdominal pressure thereby decreasing the load on the lumbar spine. Breathable elastic panels enhance patient compliance.

Available in Regular, Low Profile, Dual Panels, or Thermoformable Pockets.

Available in sizes ranging from Small to 3X-Large.

## Sacroiliac (SI) Belt



Made from an inner latex-free elastic belt with strategically placed Skinloc™ pads along with an outer nylon belt. The hi-density foam pad can be used to provide additional pressure on the sacrum.

Comfortable enough to be worn directly next to the skin.

Available in sizes Small to 2X-Large.

## EV-300 Semi-Rigid Lumbar Back Support



Lumbosacral orthosis fabricated from dense elastic material providing adjustable compression for greater lumbar containment. Incorporating a semi-rigid posterior with four malleable steel stays, threadless seam technology and a removable lumbar cushion. Also featuring a specifically designed fastening system for ease of fitting for patients with limited dexterity.

Sizes 2X-Small to 2X-Large

## Lumbitron Elite Duo Semi-Rigid Back Support



Orthosis made of breathable multi-band material, which provides uniform compression due to its threadless seam technique, avoiding annoying friction and loss of compression on the patient's skin. Reinforced in the lower back region with four preformed steel stays, which are flexible and adaptable. The front part includes 2 flexible stays for improved consistency and abdominal restraint.

Available in regular version in sizes from Small to 3X-Large as well as a short version in sizes from X-Small to 3X-Large.

## Shoulder Immobilizer



Soft padded shoulder strap provides exceptional comfort while the soft padded waist strap restricts abduction and external rotation of the shoulder. The shoulder and waist strap are now made from new material that provides greater durability as well as resistance to stretching. The new material used in the envelope is also more durable, yet very comfortable and breathes well. Bilateral design fits right or left shoulder.

Available in sizes X-Small to X-Large as well as a Universal size.

## POP Padded Arm Sling



Designed to provide an optimum level of comfort and shoulder stabilization. Features a soft, padded shoulder strap, thumb loop, stretch-resistant, lightweight, breathable material. Latex-free

Available in sizes from 2X-Small to X-Large as well as a Universal size.

## E-240 Strong Shoulder Support



Made up of two cotton elastic straps as well as independent rings for each side of the shoulders. The straps exert a gradual pull on the shoulders aiding in the correction of Kyphotic posture.

Available in sizes X-Small to 2X-Large.

## TGO330 Functional Elastic Shoulder Support



The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. In addition, it optimizes blood circulation and facilitates body temperature regulation.

Available in sizes from X-Small to X-Large

## EpiGel™ Elbow Support



A comfortable, low profile design that aids in the treatment of tennis elbow (lateral epicondylitis) or golfer's elbow (medial epicondylitis). The EpiGel elbow support provides relief with two adjustable gel pads that compress the extensor and flexor muscles without creating a tourniquet around the forearm. The soft gel pads also dampen vibration. The patented strap retention system allows for easy application.

Comes in one Universal size.

## Epitec Fix Epicondylitis Armband



A molded and shaped thermoplastic cuff lined with silicone beads adheres to the skin without applying undo pressure. The brace is made to be easily adjustable and extremely comfortable to wear.

Comes in one Universal size.

## Cubital Tunnel Brace



The Cubital Tunnel Brace is a low-profile night splint that comfortably holds the elbow in a neutral position for the treatment of cubital tunnel syndrome. Ventral design, with no exposed rigid materials, enhances patient compliance during sleep. Adjustable aluminum stay allows the elbow to be positioned at the optimum angle.

Available in sizes Small/Medium and Large/X-Large

## TGO340 TheraGo Functional Elastic Elbow Support



The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. In addition, it optimizes blood circulation and facilitates body temperature regulation.

Available in sizes from X-Small to X-Large

## OS6230 Elastic Elbow Support with Gel Pads



Controlled compression provides pain relief. The viscoelastic gel pads redistribute the compression from the prominent epicondyles to the surrounding soft tissues and provide an intermittent massage.

Sizes range from Small to Large.

## ACE401 Actius Elastic Elbow Support



Made from breathable elastic, it is lined with a highly resistant soft and stretchable fabric that provides greater comfort.

Indications include: Lateral and medial epicondylitis; Tendinitis; Mild bruising; inflammatory processes; Support and protection during sporting and work-related activities.

Available in sizes Small to Large

## 8" Wrist Lacer™ II - Wrist Support



The new 8" Wrist Lacer™ II provides equal tension across all four laces to easily generate uniform compression along the wrist support. The polypropylene felt liner and perforated suede allow moisture to evaporate away from the skin and the blue palmar pad offers enhanced support by matching the contour of the palm. The dorsal and palmar aluminum stays are removable to allow bending by a medical professional to create a custom fit.

Sizes from Pediatric to 2X-Large

## 10½" Forearm Lacer II



The new 10½" Forearm Lacer II provides equal tension across all four laces to easily generate uniform compression along the wrist support. The polypropylene felt liner and perforated suede allow moisture to evaporate away from the skin and the blue palmar pad offers enhanced support by matching the contour of the palm. The dorsal and palmar aluminum stays are removable to allow bending by a medical professional to create a custom fit.

Sizes from X-Small to X-Large

## Ryno Lacer® II Short - Wrist and Thumb Support



The new Ryno Lacer® II Short provides equal tension across all laces offering supreme comfort, ease of application and superior immobilization of the wrist and thumb. The tension can be easily adjusted for those patients with disproportionate forearm and wrist circumferences. The dorsal, palmar and radial aluminum stays are removable to allow bending by a medical professional to create a custom fit.

Available in sizes X-Small to X-Large

## Ryno Lacer® II Long - Wrist and Thumb Support



The new Ryno Lacer® II Long provides equal tension across all laces offering supreme comfort, ease of application and superior immobilization of the wrist and thumb. The tension can be easily adjusted for patients with disproportionate forearm and wrist circumferences. The dorsal, palmar and radial aluminum stays are removable to allow bending by a medical professional to create a custom fit.

Available in sizes X-Small to X-Large

## Tripod™ 8" Wrist Lacer



In addition to flexion-extension control, the Tripod Wrist Lacer resists radial-ulnar deviation. Ulnar, palmar and dorsal aluminum stays work in concert to restrict radial-ulnar deviation and minimize flexion-extension of the wrist. The polypropylene felt liner and perforated suede allow moisture to evaporate away from the skin and enhance patient compliance. The adjustable thumb strap permits unrestricted motion of the fingers.

Available in Left and Right sizes X-Small to X-Large as well as a Universal size.

## Tripod™ Lacer II - 8" Wrist Support



The Tripod Lacer II wrist support offers enhanced immobilization and ease of application. In addition to flexion-extension control, the Tripod Lacer II resists radial-ulnar deviation and its patented design provides equal tension across all laces.

Ulnar, palmar and dorsal aluminum stays work in concert to restrict radial-ulnar deviation and minimize flexion-extension of the wrist. Thumb strap has additional hook/loop interface for added security. Ulnar and palmar stays may be removed if desired for a progressive treatment protocol.

Available in Left and Right sizes X-Small to X-Large

## Compressor Wrist Support



The compressor strap on this wrist support provides compression along the carpal region. It works in conjunction with the plastic stays to create a functional wrist support that limits the range of flexion and extension of the wrist. The patented strap retention system on the end of the compressor strap allows for ease of application.

Left and Right sizes Small to X-Large.

## Motion Manager® - Wrist Support



This neoprene wrist support is an excellent choice for providing functional motion. The contoured plastic palmar stay and the multi-density dorsal stays offer adjustability to provide the optimum level of flexion and extension control. A low profile design allows the range of motion needed to return to work safely. Bilateral design fits left or right hand.

Sizes X-Small to X-Large

## OS6260 - Elastic Wrist Support



Stabilizes, supports, and provides compression to the wrist. The elastic strap can be loosened or tightened to achieve desired comfortable fit.

Bilateral design.

Available in sizes Small to Large

## OS6261 - Adjustable Wrist



Provides comfortable compression and support. Bilateral design fits both left and right sides.

One Universal size for all wrists.

## TGO350 TheraGo Functional Elastic Wrist Support



The multi-directional compression provided prevents liquid build-up, reduces stress on the soft tissues and stimulates proprioceptors, thereby enhancing articular function and muscle co-ordination. In addition, it optimizes blood circulation and facilitates body temperature regulation.

Available in sizes from X-Small to X-Large

## Allard S.O.T. Resting Hand Orthosis



S.O.T (Smart Orthotic Treatment) is a resting orthosis for the spastic and paretic hand, when the goal is to maintain or increase the mobility of the hand. Now available in pediatric to adult sizes! The orthosis is smooth and lightweight. The S.O.T has an aluminum core which allows adjustment to the desired position.

Available in Left and Right in sizes from XXX-Small to Large

## Restorative Hand



Heat moldable to provide mild stretch on the joints then remold to gradually correct joint to normal alignment. Provides prolonged low load passive stretch to re-lengthen shortened tissue. The splint is lined with Outlast temperature regulating material to not cause unwanted heat.

Color-coded straps aid in proper strap attachment.

Also available in a modified version with a dorsal attachment.

Available in Pediatric sizes from Infant to Youth as well as adult sizes from X-Small to X-Large.

## Neuroflex Restorative Flex Hand



Treats from the tightest closed palm, to moderate flexion, due to traumatic brain injury, stroke, spinal cord injury, cerebral palsy, multiple sclerosis, or Parkinson's disease.

The graduated sizes of cones gradually open hand/fingers/thumb and efficiently works to relax joint tension, even from the wrist. Reverse cones are also available to order.

Adult universal sizing in both left and right., in addition to sized braces for children.

## M760 - Manutec Fix Immobilizing Wrist w/ Palm Splint



Wrist splint manufactured from breathable velour, made using the latest generation materials (triple layer: foam, velour, sponge), with a triple immobilization system comprising a malleable aluminum palm splint and medial-lateral anatomically shaped strips for adaptation.

Indicated for wrist immobilization; post-trauma; post-surgery; painful or swollen sequelae; instability; rehabilitation.

The M760 is available sizes: Small/Medium and Large/X-Large

## SELECTION® Ulnar Deviation



This hand splint is exclusively designed and produced using our unique three layer Selection fabric. With its low profile design, it meets the contours of the hand, without impeding function. The fabric can be trimmed without fraying, and therefore allows further customization if required. The easy fastening surface allows the hook to be placed anywhere on the splint to optimize the directional pull needed. The orthosis is designed so that the material extends far to the ulnar side at 5. MCP to support the best possible direction of pull and distribute the pressure.

Available in left or right in sizes X-Small to Large

## Vission Finger Splint



Finger splint for immobilization of fingers after injury or surgery in order to prevent contractures in the finger joints. Can also be used for conservative finger therapy as well as Dupuytren syndrome and conditions requiring temporary immobilization of the fingers.

Recommended applications: Dupuytren syndrome, trauma after injury or surgery, conservative finger therapy, temporary immobilization of the fingers.

Bilateral and available in sizes from Small to Large.

## Vission Multi Finger Splint



This orthosis is intended for immobilization of fingers after injury or surgery and is aimed at preventing contractures in the finger joints. Can also be used for conservative finger therapy as well as treatment of Dupuytren syndrome. Allows temporary immobilization of your fingers. The orthosis can also be used in conjunction with treatment of fractures in the metacarpal bones 111-V and its associated Falanger.

Recommended applications: Immobilization, Dupuytren Syndrome, and fractures.

Available in left or right in sizes from Small to Large.

## Boxer Splint



The Boxer Splint is easy to adjust and offers superior comfort and support to enhance compliance. It is indicated for fourth or fifth metacarpal fractures. The malleable palmar and dorsal aluminum stays are easily bent by a medical professional to create a custom fit. The molded plastic clamshell design protects and supports the fourth and fifth metacarpals without restricting the flexion of the MP joints.

Sizes Small to X-Large

## POP Arthritis Compression Glove



The POP Arthritis Compression Gloves help to relieve the aches and pains, as well as the stiffness associated with arthritis of the hands and fingers. Our form-fitting gloves provide mild compression which helps to control the swelling and pain of the joints. The thermal lining holds in your body heat and helps to increase circulation and provide long-lasting relief. The textured palm provides a non-slip grip and open fingertips are convenient for everyday activities, allowing you to touch and feel without annoying barriers. All our gloves are latex-free. Sold by the pair.

Available in sizes XX-Small to XX-Large.



## M780 Manutec Fix Rizart Plus CMC Thumb Brace



Highly resistant and very soft material. The brace is heat-moldable, flexible and comfortable, and perfectly adapts to the shape of the hand.

The design with double injection of materials and the rounded edges guarantee maximum comfort of use while the fastening system, secured by a metal core, enables the degree of stability of the trapeziometacarpal joint to be adjusted as required.

Available in left or right in sizes Small to Large.

## CMCcare Thumb Brace



Contours around the first metacarpal to stabilize the CMC joint for pain-free motion.

The CMCcare restores proper alignment of the thumb for improved function.

Available in Left and Right sizes from Small to Large.

## TeePee™ Thumb Protector



Three aluminum stays comfortably support the MCP joint. Each stay can be shaped and contoured to provide a custom fit. Fits either left or right thumb.

Bilateral design available in sizes X-Small to X-Large.

## M770 Breathable Thumb Immobilizing Splint



Made of breathable velour. Consists of a malleable aluminum splint which can be molded to provide a custom fit. Latex Free. Bilateral design can be adapted for both the left and right hand.

Available in Small and Large

## Universal Suede Thumb Support



Easy to apply, comfortable and low profile splint. Anatomically contoured aluminum stay and soft elastic closure straps comfortably immobilize the MP & CMC Joints of the thumb.

Bilateral design comes in a short and long version in one Universal.

## CMC-X Thumb Stabilizer



The CMC-X™ Thumb Stabilizer provides your patients with pain relief from carpometacarpal (CMC) joint arthritis while allowing them to continue the activities they enjoy.

Available left or right in sizes Small to X-Large.



Professional Orthopedic Products

Thank you for choosing Professional Orthopedic Products for your orthopedic, neurological, and sports medicine product needs. We appreciate and value your business and have assembled the following to make your ordering experience with Professional Orthopedic Products as simple and clear as possible.

### **OUR MISSION STATEMENT**

At Professional Orthopedic Products, we are dedicated to providing high quality orthopedic, neurological, and sports medicine products while maintaining high standards of customer service.

### **DEDICATED CUSTOMER SERVICE**

Telephone:	(514)-484-4441	Toll Free:	1-866-334-4441
Fax:	(514)-484-4431	Toll Free Fax:	1-888-484-4431
E-mail:	csr@poportho.com		

All orders should be faxed or e-mailed to ensure accuracy. For your convenience, we can supply you with an order form upon request. In most cases, orders received before 3:00 PM EST will be shipped out the same day.

### **WARRANTY**

All products sold by Professional Orthopedic Products have a limited warranty against defects in materials or workmanship for a period of six months from date of purchase. Before returning a product you believe is defective, you must contact Professional Orthopedic Products for a Return Authorization number. After receipt and inspection of the product, should the product prove to be defective, we will repair or replace it at no charge. This warranty is limited to materials and workmanship only, and does not include failure due to wear, incorrect laundering, or improper fitting. Please note that all the returns **MUST** be accompanied by a dated proof of purchase to be eligible for a credit or exchange.

### **RETURN POLICY**

Before returning or exchanging any product for reasons other than a defect, you must contact Professional Orthopedic Products for a Return Authorization number. This will insure prompt and accurate credit or replacement. Returns and exchanges will only be accepted within sixty days from the date of purchase and will be subject to a 25% restocking fee. Custom products, discontinued products, products that have been changed in material or design or products that have been worn may not be returned or exchanged for any reason other than POP's error. Please note that all the returns for reasons other than defects must be accompanied by a dated proof of purchase and must be returned with all the original packaging intact to be eligible for a credit or exchange.

### **CREDITS**

Credits will be mailed to you as soon as they are issued. For your convenience, we leave it up to you to inform us as to when and to which invoice you would like the credit applied.

In closing, we look forward to serving you and to a long and prosperous relationship.

**Professional Orthopedic Products**

# Professional Orthopedic Products



Professional Orthopedic Products carries the following brands:




Elite Orthopaedics



Professional Orthopedic Products

*Providing Innovative Products  
of Unsurpassed Quality*



Telephone 1-866-334-4441 • Fax 1-888-484-4431 • Email [csr@poportho.com](mailto:csr@poportho.com)

[www.poportho.com](http://www.poportho.com)