



EP-21G EPITEC FIX EPICONDYLITIS ELBOW BAND

Features:

The exterior is made with highly resistant cushioned materials while the interior features a molded and shaped thermoplastic body with silicone points that adhere to the skin by applying pressure. Easily adjustable with the closing band and the elastic band.

Function:

Unloading, compression and protection.

Bilateral design fits both left and right elbow

ITEM	SIZE
EP-21G	UNIVERSAL

Indications:

- Prevention of sporting and work-related injuries
- Lateral Epicondylitis (tennis elbow)
- Medial Epicondylitis (golfer's elbow)



With the EPITEC FIX straps, pressure is applied through a viscoelastic cushion with special adjustment points. Correctly positioned, these reduce the force of traction which results from the insertion of the epicondyle musculature tendon; achieving this through the orientation of the muscular fibres.