

# EpiGel<sup>®</sup> with Stays Elbow Support



The new EpiGel with Stays elbow support has enhanced focal compression to aid in the treatment of tennis elbow (lateral epicondylitis) or golfer's elbow (medial epicondylitis). The EpiGel with Stays provides relief with two adjustable gel/stay pads that compress the extensor and flexor muscles without restricting circulation of the forearm. The universal design fits a broad size range of patients.



### Two gel/stay pads

The polymer gel and aluminum stay pads are connected with hook and loop fastener and work in concert to compress extensor muscles to relieve lateral epicondylitis and flexor muscles to relieve medial epicondylitis. The gel/stay pads also elevate the strap off the skin where there is venous flow return so there is no tourniquet around the forearm. The contour of the aluminum stays may be adjusted, and the gel/stay pads may be repositioned to provide a custom fit. The soft gel pads also dampen vibration.

### Unique strap retention system

Allows for easy application.

### Universal design

Fits right or left arm. One size fits most.

### Low profile design

Less likely to catch on any object.

Item #	Size	Elbow Circ.
223440	univ	8" - 14"

\*measure forearm just below elbow

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060  
p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com

EpiGel and Med Spec are trademarks of Medical Specialties Inc.

© 2024 Medical Specialties Inc.

